

SIGNA PACKAGE




REGIUS XII



IMAGING STATION
IMAGEPILOT
Touch



DRYPRO 2




Giving Shape to Ideas



Konica Minolta has been a CR technology leader since 1996-consistently designing, developing, and delivering innovative solutions. The REGIUS ΣII represents the sum of this experience. REGIUS ΣII delivers a compact and affordable imaging solution that maintains the superior image quality and reliability expected from our REGIUS family of products. Combined with the redesigned ImagePilot software, the REGIUS ΣII provides simple operation and smooth workflow to a wide range of medical practitioners.


High Throughput of 60plates/hour(maximum)

REGIUS **SII**



From patient registration to long term archive





I. Register

Manual or automated patient registration from anywhere in the office. Select the patient from the worklist or input patient information. No need to specify body part or exam type.

II. Digitize

Just insert an exposed cassette into the REGIUSΣII CR reader. The ImagePilot then automatically processes the image for optimized diagnostic viewing.

ImagePilot

Quick Preview

The fully automated image processing produces consistent optimized images that display within 45 seconds



Print composition


ImagePilot supports a print composer mode. Use this mode from any client station to print multiple images on one sheet of film or paper.



▲Print Composer

Easy recovery

If the plate jam occurs, it can be removed easily by opening the top cover or the backside cover.



Multifunctional Viewing

ImagePilot supports, over 100 features and functions.

·Annotation

Line/Arrow/Rectangle/Text/ Shutter/Marker, etc

·Scaling

Cardio Index/Distance/ Angle/Rectangular CT/ Pelvis Analysis, etc

·Image Tool

Window Level/Zoom/Pan Magnifier/ REGIUS Processing, etc

Report

Other Modality Image






▲Report

■Web Viewing

At your home and trip destination, you can see the images of image Pilot via VPN.



Konica Minolta provides complete workflow





II. View

Images can be viewed at up to 20 ImagePilot client workstations throughout the office. Each client supports all the same ImagePilot functions - the same GUI at every seat.


IV. Print

The DRYPRO SIGMA outputs the films up to 45 sheets (14x7")per after the image has been printed from ImagePilot DRPRO SIGMA has four types of films available, which can be transposed based on user needs.

Drypro Sigma

Compact Design

DRYPRO SIGMA is a desktop DRY Laser Imager. It features a compact footprint as small as 65x60cm. It is designed to seamlessly integrate into various environments when combined with REGIUS $\Sigma_{\rm II}$



Four film sizes and easy loading process

Users can use four types of films for different purposes. Open the front cover, retrieve the film inside, and load new one. Replacing the film can be done with this simple process




User – Friendly

DRYPRO SIGMA is easy to operate, requiring no hands – on training upon installation and offering intuitive usability.

DICOM Printer

DRYPRO SIGMA can connect easily to other CR and DR consoles. CT, and MR, in addition to ImagePilot. Non-DICOM images can be scanned easily using imagePilot. Non-DICOM images can be scanned easily using ImagedPilot and can also be printed using DRYPRO SIGMA.



COMPACT

REGIUS ΣII

Konica Minolta's Ultra-Compact CR

The REGIUSΣI is Konica MinoIta's smallest and lightest CR. The desktop CR has a footprint of only 0.31m² and weighs a mere 28kg.

0.11m³ & 28kg



Environmentally Friendly

The REGIUS ΣI was designed with the environment in mind. Energy consumption has been reduced to 100VA - about the

same as a household light bulb. That's one tenth the power of the smallest film processor.



HIGH QUALITY


High Resolution

The REGIUS Σ I supports pixel sizes of 87.5 μ m and 175 μ m for all cassette sizes, achieving a maximum resolution of 4020 \times 4892 pixels (14 \times 17 at 87.5 μ m).



Integral Processing

Integral Processing - the result of Konica Minolta's vast imaging experience - is a robust algorithm that produces the optimum diagnostic image automatically. This patented technique greatly simplifies x-ray image acquisition while delivering incredibly consistent image quality and virtually eliminating retakes.



SMART CONSOLE

ImagePilot

All-in-One!

The ImagePilot software combines CR image acquisition and PACS functionality. It is a complete workflow solution for clinics, imaging centers, and small hospitals. From one GUI, users can acquire, review, distribute, and archive patient images.



Easy Control

ImagePilot supports a multi-touch display for the easiest operation. The more than 100 functions and features are literally at the tips of your fingers.



REGIUS Σπ Specifications

Cassette size

14"×17" / 14"×14" / 11"×14" / 10"×12" / 8"×10" / 18×24cm / 24×30cm / 15×30cm

Sampling Pitch

2Types: 87.5μm / 175μm

■ Maximum Resolution

4,020×4,892 (14"×17", 87.5µm)

■ Digital Gradation Level

4096 levels (12bit)

Processing Capability

Up to 60 plates / hour (14"×17" / 175 μ m) Up to 45 plates / hour (14"×17" / 175 μ m)

■ Slots

Insertion / ejection slot × 1

Outer Dimensions/Weight

W510×D610×H355mm/Approximately 28kg

Power Supply Conditions

100VA max during operation: 20VA when not in use.

Operating Conditions

Temperature: 10-30℃

Humidity : 15-80%RH (No Condensation)

Readable Cassette

REGIUS CASSETTE RC-300

14"×17", 14"×14", 11"×14", 10"×12", 8"×10", 24×30cm, 18×24cm, 15×30cm

Readable Plate

REGIUS PLATE FP-1S

14"×17", 14"×14", 11"×14", 10"×12", 8"×10", 24x30cm, 18×24cm, 15×30cm

DRYPRO SIGMA Specifications

Laser Source

Semiconductor Laser

Film Size

 $14"\times17"(35\times43cm), 11"\times14"(28\times35cm), 10"\times12"(25\times30cm),$ and $8"\times10"(20\times25cm)$ selectable

Film

Dry Imaging recording film SD-S

Image Format

1,2,4,6,8,9,12,15,16,20,24,30,35

■ Image Memory

2GB

Pixel Size

78 u m

Processing capability

Up to 45 films per hour (14x17in. (35x43cm))

Input interface

Ethernet 10/100/1000 Base-T

- Integrated DICOM interface supports printing from DICOM print modalities
- *Connect non-DICOM modalities using PrintLink IV,ID

Operating condition

15 - 33°C (59-91F) 20 - 80% RH, non-condensing

Power

100-120V \pm 10% (50/60HZ) 9A 200-240V \pm 10% (50/60HZ) 5A


Dimensions / Weight


W597mm X D649mm X H463mm / 54kg

Footprint

0.39m²


REGIUSΣII Outer Dimensions






Unit:mm 1/10 scaled

DRYPRO SIGMA Outer Dimensions





Specification are subject to change without prior notice

*Specification are subject to change without prior notio





KONICA MINOLTA HEALTHCARE INDIA PVT. LTD.

Office No.515, 5th Floor, C-wing, 215-Atrium Andheri (east), Mumbai 400 059. INDIA. Tel.: +91 - 22 - 6196969 Fax: +91 - 22 - 61916996

email: sales@mi.konicaminolta.in www.konicaminolta.in/healthcare