

CE IVDR

azooka



# AZUL PLASMID DNA EXTRACTION KIT

DNA IN 40 MINS | GOOD YIELDS FOR USE IN PCR/SEQUENCING

## PRODUCT BROCHURE



Cat No-DE104

ISO 13485 CERTIFIED

## PRODUCT DESCRIPTION

AZUL Plasmid DNA Extraction Kit is an easy and efficient system for the isolation of plasmid DNA from bacterial cells. This kit uses a silica-based spin column technology for isolating total DNA from biological samples. The eluted plasmid DNA is suitable for Cloning, transfection, transformation, restriction endonuclease digestion in vitro transcription/translation, qPCR, and Next-Generation sequencing.

## KIT COMPONENTS

| Components         | For 50 preps    | For 25 preps    |
|--------------------|-----------------|-----------------|
| P1 Buffer          | 15mL            | 8mL             |
| P2 Buffer          | 15mL            | 8mL             |
| P3 Buffer          | 20mL            | 10mL            |
| Binding buffer(BB) | 25mL            | 13mL            |
| Wash Buffer (WB)   | 50mL            | 25mL            |
| Elution Buffer(EB) | 4mL             | 2mL             |
| Spin Column        | 50 (Pouch pack) | 25 (Pouch pack) |

## SPECIFICATIONS


|                   |   |
|-------------------|---|
| Format            | Spin column   |
| Sample type       | Bacterial cells   |
| Equipment         | Microcentrifuge   |
| Processing time   | <40 mins  |
| Processing volume | 1mL   |
| Type              | Plasmid DNA   |
| Sample storage    | Eluted DNA should be stored at $\leq -20^{\circ}\text{C}$ |
| Yield             | 80 ng - 700ng/ $\mu\text{L}$                              |
| Purity            | $A_{260}/A_{280} \geq 1.8$                                |
| Kit Storage       | Room Temperature  |
| Kit Validity      | Viable for 1 year if stored at appropriate conditions     |

**NOTE:** Check the Buffers for any salt precipitation before every use. Re-dissolve any precipitate by warming the solution to  $37^{\circ}\text{C}$ , then cool it back to room temperature before use.

## DNA EXTRACTION PROTOCOL

1. In a 1.5mL microfuge tube, transfer up to 1mL of the bacterial culture and centrifuge at 15,000 rpm for 2 mins to pellet the cells.
2. To the pellet, add 200µL-300µL of P1 Buffer. Mix briefly by vortexing for 30 seconds.
3. Add 200µL-300µL of P2 Buffer to the same tube and mix by inverting gently until the solution turns transparent. Do not vortex. Incubate at RT for 5 min.
4. Add 350µL of P3 Buffer and mix by inverting 3-4 times. Centrifuge at 15,000 rpm for 10 min.
5. Carefully transfer the clear supernatant to the new 1.5 mL microfuge tube. 500µL of Binding Buffer is added to this tube and mixed slowly by inverting the tube 5 times. Incubate at -20°C for 15 min.
6. Transfer up to 800µL lysate to the spin column inserted in a collection tube. Centrifuge at 13,000 rpm for 3 min. Discard the flow-through and place the purification column back into the collection tube.  
Repeat this step until the entire lysate has been transferred into the column and centrifuged.
7. Add 500µL of Wash Buffer (WB) to the column and centrifuge at 15,000 rpm for 1 min. Discard the flow-through and place the purification column back into the collection tube.
8. Repeat step 7 once again to remove salts and impurities completely.
9. Spin the column at 15,000 rpm for 1 min to dry the column. Place the tube at RT for 2 mins.
10. Transfer the purification column to a clean, sterile microfuge tube and add 50µL of Elution Buffer or DNase/RNase-free water to the centre of the column. Centrifuge the column for 15,000 rpm for 1 min.
11. Discard the purification column and store the eluted Plasmid DNA at -20°C or -80°C until use.

## FLOW DIAGRAM OF DNA EXTRACTION PROTOCOL



## TROUBLESHOOTING

| PROBLEM                   | POSSIBLE CAUSES  | SUGGESTED SOLUTIONS  |
|---------------------------|--|--|
| Low DNA Yield             | <b>Too much culture used:</b> Incomplete lysis and neutralization are two of the most common causes of failed plasmid preps, and both are caused by too much culture being used. | Use less input material or increase the volume of the Lysis Buffer and neutralization buffer accordingly.  |
| Low DNA Purity(A260/A280) | Improper Sample Handling causes Ethanol or Salt contamination.   | Make sure lysate and wash buffers have passed entirely through the matrix of the column. This may require centrifuging at a higher speed or longer time. |
| Genomic DNA Contamination | Improper Sample Handling: Sample was vortexed or handled too roughly.  | Slightly invert the tubes for mixing.  |
|                           | Genomic DNA contamination is usually caused by excessive mechanical shearing during the lysis and neutralization steps.  | Avoid vortexing after adding lysis and neutralization buffer.  |
|                           | Prolonged or Incomplete lysis or neutralization may contribute to genomic DNA contamination.   | Follow protocol to avoid such mistakes.  |
|                           | <b>Overgrown or Old culture:</b> May contain more genomic DNA. contamination than fresh cultures.  | Use fresh culture for optimal performance.   |
| RNA Contamination         | Too much culture used.   | Reduce the volume of culture being processed.  |
|                           | Incorrect column washing.  | Use the wash volume mentioned in the protocol.   |

## ORDERING INFO

| CATALOG NO | PRODUCT  | PREP        |
|------------|--|-------------|
| DE101      | AZUL Tissue DNA Extraction Kit                 | 25/50 preps |
| DE102      | AZUL Animal Cell Culture DNA Extraction Kit    | 25/50 preps |
| DE103      | AZUL Bacterial DNA Extraction Kit              | 25/50 preps |
| DE104      | AZUL Plasmid DNA Extraction Kit                | 25/50 preps |
| DE105      | AZUL Plant DNA Extraction Kit                  | 25/50 preps |
| DE106      | AZUL Soil DNA Extraction Kit                   | 25/50 preps |
| DE107      | AZUL Blood DNA Extraction Kit                  | 25/50 preps |
| DE108      | AZUL Cell-free DNA Extraction Kit              | 25/50 preps |
| DE109      | AZUL DNA Extraction Kit- Difficult samples     | 25/50 preps |
| DE110      | AZUL Saliva DNA Extraction Kit                 | 25/50 preps |
| DE111      | AZUL Stool DNA Extraction Kit                  | 25/50 preps |
| DE112      | Quick AZUL Bacterial/Fungal DNA Extraction Kit | 25/50 preps |
| DE113      | AZUL Microbiome DNA Extraction Kit             | 25/50 preps |
| DE114      | AZUL Gel DNA Extraction Kit                    | 25/50 preps |
| DE115      | AZUL FFPE DNA Extraction Kit                   | 25/50 preps |
| DE116      | AZUL Chloroplast DNA Extraction Kit            | 25/50 preps |
| DE117      | AZUL Mitochondrial DNA Extraction Kit          | 25/50 preps |
| DE118      | AZUL Pollen DNA Extraction Kit                 | 25/50 preps |
| DE119      | AZUL Fungal DNA Extraction Kit                 | 25/50 preps |
| DE120      | AZUL Sperm DNA Extraction Kit                  | 25/50 preps |
| DE121      | AZUL Skin DNA Extraction Kit                   | 25/50 preps |

## FEEDBACK

### How did this kit perform?

Did AZUL Extraction Kit fulfill expectations required for your research?

Let us know by filling out the feedback form [here](#)

Or scan the QR code:



## CONTACT US



hello@azooka.life



+91 8088747968



www.azooka.life



AZOOKALIFE



AZOOKA.LABS



AZOOKALIFESCIENCES



AZOOKALIFE



### RESEARCH CENTRE:

Society for Innovation and Development, Indian Institute of Science, Malleshwaram, Bengaluru, Karnataka, India- 560055

### MANUFACTURED AT:

# 1A, Kushal Garden Arcade, 'C' Block, 5th Floor, Peenya Industrial Area, 2nd Phase, Bengaluru, Karnataka, India- 560058