

INTENDED USE

AZUL Pollen DNA Extraction Kit is an easy and efficient system for the isolation of pollen DNA from samples like honey.

SUMMARY AND EXPLANATION

This kit uses a silica-based spin column technology for isolating DNA from biological samples, thereby eliminating toxic phenol-chloroform extractions. The eluted DNA is suitable for all sensitive downstream applications such as qPCR and Next-Generation sequencing.

PRODUCT FEATURES

- Rapid purification of high-quality, ready-to-use DNA.
- No organic extraction or alcohol precipitation.
- Consistent and high yields.
- Complete removal of contaminants and inhibitors for reliable results.
- Kit formats for low- to high-throughput – options for automation of all kits.

PRECAUTIONS

- Avoid all skin contact with reagents in this kit. In case of contact, wash thoroughly with water.
- AZUL Pollen DNA Extraction Kit is intended for use as supplied. Do not dilute or add other components to the AZUL Pollen DNA Extraction Kit.

DIRECTIONS FOR USE

- 1.Take 3 mL of Honey ( which might contain pollen cells) dissolve in 1 ml sterile water and incubate at 65°C for 30 min.
- 2.Centrifuge at 5000 rpm for 10 min. Discard the supernatant, and air dry the pellet for 5 mins at room temperature.
- 3.Add 700 µl extraction buffer 1, 4 glass beads and vortex well for 1-2 mins.
- 4.Add 100 µl of lysis buffer 1 and 10 µl proteinase K, mix by gentle inversion, and incubate at 56°C for 1 hour.
- 5.To this add 500 µl extraction buffer 2, 10 µl proteinase K, and 50 µl lysis buffer 2, mix well and incubate at 65°C for 4 hours - overnight.
- 6.Centrifuge at 10,000 rpm for 10 mins. Transfer the clear supernatant to a 2 ml Eppendorf tube, further add 500 µl of binding buffer and mix well.
- 7.Transfer the lysate to a clean spin column. Centrifuge the spin column at 13,000 rpm for 2 min at RT.
- 8.Discard the flow-through and place the purification column back into the collection tube. Repeat this step until the entire lysate has been transferred into the column and centrifuged.
- 9.Wash the spin column with 500 µL Wash Buffer (WB) at 13,000 rpm for 1 min and discard the flow through. Repeat this step again.
- 10.Keep the purification column in a clean, sterile 1.5 mL microfuge tube and add 30 µL of Elution Buffer or DNase/RNase-free water to the center of the column.
- 11.Centrifuge the column for 13,000 rpm for 2 min.
- 12.Discard the purification column and store the eluted DNA at -20°C or -80°C until use.

KIT COMPONENTS

| Components          | For 50 preps    | For 25 preps    |
|---------------------|-----------------|-----------------|
| Extraction Buffer 1 | 36mL            | 18mL            |
| Lysis Buffer 1(LB1) | 6mL             | 3mL             |
| Extraction Buffer 2 | 26mL            | 13mL            |
| Lysis Buffer 2(LB2) | 3mL             | 1.5mL           |
| Proteinase K        | 1mL             | 0.5mL           |
| Glass Beads         | 200             | 125             |
| Binding buffer(BB)  | 30mL            | 15mL            |
| Wash Buffer (WB)    | 60mL            | 30mL            |
| Elution Buffer(EB)  | 4mL             | 2mL             |
| Spin Column         | 50 (Pouch pack) | 25 (Pouch pack) |

CAUTION

- Check the buffers for any salt precipitation before every use.
- Re-dissolve any precipitate by warming the solution to 37°C, then cool it back to room temperature before use.
- During operation, always wear a lab coat, disposable gloves, protective goggles and mask.

KIT STORAGE AND STABILITY

- Store the kit at room temperature.
- Viable for 1 year if stored at appropriate conditions.

ORDERING INFORMATION

Please call us at +91 8088747968 or mail at hello@azooka.life for any queries or assistance. Additional information can be found online at www.azooka.life