

Features:

- Color OLED display
- Multi-direction display
- Multi-mode display
- PR waveform and bargraph display
- Accurate measurement of SpO₂ value and PR value
- Power off automatically when no signal
- Audible and visible alarm function
- Low battery voltage indication
- Low power consumption
- Widely used in hospital, home, community medical treatment, sports healthcare etc.

Specifications

Model			BP-10M		
Display			Color OLED		
Display Direction			4		
Waveform			YES		
Audio Alarm			YES		
SpO ₂	Measurement Range		35%-100%		
	Accuracy		70%-100%, ±2%; 35%-69% no definition		
	Resolution		1%		
PR	Measurement Range		25-250 bpm		
	Accuracy		± 2 bpm		
	Resolution		1 bpm		
Low Perfusion			≤0.3‰		
Alarm Indication	SpO ₂		< 90%		
	PR	Upper Limit Lower Limit	5-250 bpm (default 130 bpm) 5-250 bpm (default 50 bpm)		
Power Supply Management			Power off automatically (no signal for more than 16 sec)		
Power Requirements			AAA size battery (1.5V)×2		
Battery Life			24 hrs		
Power Consumption			30 mA		
Dimension			64×38×37mm (L×W×H)		
Net Weight Including Battery			57 g		

Shenzhen Biocare Bio-Medical Equipment Co., Ltd.


Headquarters: 7/F, Block A, Shenzhen Mingyou Industrial Products Exhibition & Procurement Center, Baoyuan Road, Xixiang Sub-district, Bao'an District, 518102 Shenzhen, P.R. China

Factory: 2/F West, 4th Block, Dayang Road South, Fuyong Sub-district, Bao'an District, 518103 Shenzhen, P.R. China

Tel: +86-755-3661 5333 **Fax:** +86-755-2796 0643

E-mail: sales@biocare.com.cn

© 2014 BIOCARE all rights reserved.
ENG-CATA-V1.0-20140928



BP-10M

Fingertip
Pulse Oximeter



www.biocare.com.cn