

Features:

- Elegant and lightweight design
- Color TFT LCD display of SpO₂, PR, Pulse bar, PI and Plethysmogram
- Real-time review SpO₂ ang PR value
- Up to 70-hour of oximeter data storage (Record data once every 1/2/4/8 seconds)
- Data record can be uploaded to PC
- Oximeter data manager software is available
- Adible & visible alarm function
- Low battery voltage indication
- Power off automatically when no signal
- More than 40 hours battery life
- AC/DC adapter as optional



Specifications

| | | |
|------------------|---|---|
| Model | | BP-12C |
| SpO ₂ | Measurement range | 35%-100% |
| | Accuracy | ≤ 3% (35%-100%) |
| | SpO ₂ low over-limit setting range | 85%~99% |
| | Default low over-limit | 85% |
| PR | Measurement range | 30 bpm~240 bpm |
| | Accuracy | ± 2% or ± 2 bpm (whichever is greater) |
| | PR over-limit setting range | 25 bpm~250 bpm |
| | Default over-limit | 120 bpm (high), 50 bpm (low) |
| Perfusion Index | Measurement range | 0~20% |
| Power supply | Battery type | 3 x LR6 (AA) alkaline batteries or Ni-MH rechargeable batteries |
| | Supply voltage | 4.5V |
| | Operating current | < 50mA |
| | Continuous working time | > 30 hours |
| Data Record | | Record data once every 1/2/4/8seconds, up to 70–hour records |
| Dimensions | | 145 mm (L) × 74 mm (W) × 29 mm (H) |
| Net Weight | | 210g (including batteries) |



www.biocare.com.cn

BP-12C
Handheld Pulse
Oximeter



Shenzhen Biocare Bio-Medical Equipment Co., Ltd.

Headquarters: 7/F, Block A, Shenzhen Mingyou Industrial Products Exhibition & Procurement Center, Baoyuan Road, Xixiang Sub-district, Bao'an District, 518102 Shenzhen, P.R. China

Factory: 2/F West, 4th Block, Dayang Road South, Fuyong Sub-district, Bao'an District, 518103 Shenzhen, P.R. China

Tel: +86-755-3661 5333 **Fax:** +86-755-2796 0643

E-mail: sales@biocare.com.cn

© 2015 BIO CARE all rights reserved.
ENG-CATA-V1.0-20150415

