



## AQUA-XL Water Analysing Kits

### Phosphate

Code : XL 501

#### Directions for use :

- 1) Take 5 ml of water sample with help of dropper in 10 ml glass jar.
- 2) Make it up to 10 ml mark with Distilled or Phosphate free raw water.
- 3) Add 10 drops of Reagent LP1, Mix Well.
- 4) Add 2 drops of Reagent LP2, Mix Well.
- 5) Wait for 5 minutes.
- 6) Place the glass test jar on colour card & match with different colours. (View from top side-vertically downwards).
- 7) Note down the ppm level of O-Phosphate as  $\text{PO}_4$  after colour comparison.



## AQUA-XL Water Analysing Kits

### Phosphate

Code : XL 511

#### Directions for use :

- 1) Take 1.2 ml of water sample in the 10 ml glass jar. Make it up to 10 ml mark with Distilled or Phosphate free raw water.
- 2) Add 10 drops of Reagent LP1. Mix Well .
- 3) Add 2 drops of Reagent LP2, Mix Well.
- 4) Wait for 5 minutes.
- 5) Place the glass test jar on colour card & match with different colours. (View from top side-vertically downwards).
- 6) Note down the ppm level of O-Phosphate as  $\text{PO}_4$  after colour comparison.

