

KALAM FLUOR

Product Information

Kalam Fluor is an eco-fluorophore protein gel stain for rapid colorimetric protein band detection. Instantly view intensely colored green fluorescent bands under UV illumination. Kalam Fluor is an equal alternative to Coomassie Brilliant Blue R-250 and is compatible with all existing gel documentation systems.




Fig. 3: Protein bands stained with Kalam Fluor

Product Features

- Detection time: 1-5 min
- Sensitivity: Matches Coomassie R-250
- Proteins can not be further processed for blotting or sequencing after staining

Ordering Information

Cat. No.	Product Name	Unit Size	Price
KF101-A	Kalam Fluor for Proteins	250 mL	1094.5
KF101-B	Kalam Fluor for Proteins	500 mL	2189